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INDEX

PRESENTATION OF MANUAL	23
IDENTIFICATION OF EQUIPMENT	24
Identification	24
Description of the Equipment	25
Indication of Equipment	25
MODULES, ACCESSORIES, OPTIONS AND MATERIALS OF CONSUMPTION	26
TECHNICAL SPECIFICATIONS	27
- Technical characteristics	27
- Security and protection system	28
- Eletromagnetic emissions	28
- Norms employed:	31
- List of pieces and circuit scheme	32
- Dimensions	32
- Symbologies of packaging	33
- Symbologies of product	33
- Contents of the accessible and inaccessible markings	33
OPERATION OF EQUIPMENT	34
PRECAUTIONS, RESTRICTIONS AND WARNINGS	35
- Transportation, storage and operation	35
- Sensitiveness to environmental conditions foreseeable in normal situations of use	35
- Precautions and warnings "during the installation" of equipment	35
- Recommendations for preserving the equipment	36
- Precautions and warnings "during the utilization" of equipment	36
- Precautions and warnings "after" the use of equipment	36
- Precautions and warnings during the "cleaning and disinfection" of equipment	37
- Precautions in case of alteration in the functioning of equipment	37
- Precautions to be adopted against foreseeable or uncommon risks, related to the deactivation and abandoning of equipment	37
CORRECTIVE AND PREVENTIVE MAINTENANCE AND PRESERVATION	37
- Preventive maintenance	37
- Corrective maintenance	37
- Additional procedures for reuse	37
- Cleaning and disinfection	38
UNFORESEEN EVENTS – SOLUTION OF PROBLEMS	39
WARRANTY OF EQUIPMENT	40
FINAL CONSIDERATIONS	40

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ATTENTION

For greater safety:

- Read and understand all the instructions contained in these Instructions for Use before installing or operating this Equipment.

Note: These Instructions for Use must be read by all the operators of this Equipment.

IDENTIFICATION OF THE EQUIPMENT

Dear Customer

Congratulations. You have made a good choice when you decided to buy a GNATUS QUALITY product comparable to the best products available in the World. This manual is a general presentation of your product and it will give you important details to help you to solve possible problems.

Please, read it and keep this with you.

Identification

Technical Name: Controller for Electric Micromotor

Trade name: Electrical Micromotor

Model: Celeritá

Brand: Gnatus

IDENTIFICATION OF THE EQUIPMENT

Description of Equipment

Electrical Micromotor extremely versatile, used by demanding professionals that want an agile equipment for aesthetic use.

Its design offers light weight and convenience, ease of movement, simplicity of handling.

- Bivolt;
- Provided with safety device with anti-heat;
- Gradual adjustment of speed from 0 to 20.000 rpm;
- Direction of adjustable rotation (left/right);
- It allows the use of many types of drills and emery board;
- Low level of noise and vibration;
- Lubrication is not necessary.


Indication of Equipment

Electrical Micromotor developed for podiatry, removes and cleans the cuticle, preparation of foot and nails, for works of wearing, polishing, silicone nail finishing, gel, porcelain and polygel.

It should be used and handled by trained person (and professional personnel regulated according to local legislation of the country) following the instructions contained in this manual.

It is the user's obligation to use only the equipment in perfect condition and protect themselves, patients and others against the hazards.

MODULES, ACCESSORIES, OPTIONS AND MATERIALS OF CONSUMPTION

- 01 - Potency Controller
- 02 - Electrical micromotor hand piece
- 03 - Rotation inversion selector
- 04 - Micromotor connection
- 05 - Switch
- 06 - Voltage selector (110/220V)
- 07 - Hand selector (H) / foot (F)
- 08 - Foot connection
- 09 - Power cable
- 10 - LED indicator ON/OFF
- 11 - Bit
- 12 - Sandpaper (not included)
- 13 - Foot control (optional)
- 14 - Support handpiece (optional)

The content of this page is for information purposes, and the device may be presented different from that illustrated. Therefore, when purchasing the product check the technical compatibility between the device, coupling and accessories.

TECHNICAL SPECIFICATIONS

Device's Technical Features and Accessories

General

Classification of Equipment

According to the norm NBR IEC 60601-1

Electric shock protection type

Class II Equipment

Degree of protection against electric shock

B type

Protection against harmful water penetration

IPX 0

Mode of Operation

Continuous with intermittent load: T-on 10min. - T-off 10min

Electrical Micromotor – Control Unit

Input Voltage:

Bivolt :127/220V

Output voltage:

DC 24V 1.35A

Frequency:

50/60 Hz

fuse:

2A

Power:

50VA

weight (g):

882

Electrical Micromotor – Handpiece

Input Voltage:

DC 24V 1.35A

Torque:

2.7 N/cm

maximum rotation (RPM):

20.000

weight (g):

162

TECHNICAL SPECIFICATIONS

Security and protection system

In order to protect the engine and the control unit from the excessive temperature caused by overload, such as bearing failure, load in excess, etc., a temperature sensor is activated in a way that when the temperature overcome a certain pre-adjusted value, the Power supply is turned off.

When this protection system is operating, turn the unit off and, after 10 seconds, turn it on again.

Eletromagnetic emissions

Recommended distances between portable and mobile RF communication equipments and the Celeritá Electrical Micromotor			
The Celeritá Electrical Micromotor is made to be used in an electromagnetic environment in which RF disturbances are controlled. The client or the user of the Celeritá Electrical Micromotor may help preventing electromagnetic interference by keeping a minimal distance between mobile and portable RF communication equipment (transmitters) and the Celeritá Electrical Micromotor, as recommended below, in accordance with the maximal voltage output of the communication equipment.			
Transmitter Maximum Output (W)	Separation distance according to transmitter frequency (M)		
	150 kHz to 80 Mhz d= 1,2√p	80 kHz to 800º Mhz d= 1,2√p	800 kHz to 2,5º GHz d= 2,3√p
0,01	0,12	0,12	0,23
0,1	0,38	0,38	0,73
1	1,2	1,2	2,3
10	3,8	3,8	7,3
100	12	12	23
For transmitters with a maximum nominal output power not listed above, the recommended d separation distance in meters (M) can be determined using an equation applicable to the frequency of the transmitter, where P is the transmitter maximum nominal output in watts (W) according to the transmitter manufacturer.			
NOTE 1 At 80 MHz and 800 MHz, is applied the separation distance for the higher frequency range.			
NOTE 2 These guidelines may not apply to all situations. The absorption and reflection from structures, objects and people affect the electromagnetic propagation.			

TECHNICAL SPECIFICATIONS

Guidelines and manufacturer's declaration - electromagnetic immunity			
The Celerit� Electrical Micromotor is made to be used in the electromagnetic environments specified below. The client or the user of the Celerit� Electrical Micromotor must be sure that it is used in such environment.			
Immunity test	ABNT test level NBR IEC 60601	Level of compliance	Electromagnetic Environment Directives
RF conducted IEC 61000-4-6 RF radiated IEC 61000-4-3	3 vrms 150 kHz up to 80 MHz 3 V/m 88 MHz up to 2,5 GHz	3 Vrms 3 V/m	It is advisable that portable and mobile RF communication equipment is not used near any part of the equipment, including cables, with a separation distance less than the one recommended, calculated from the equation applicable to the frequency of the transmitter. Recommended separation distance: $d = 1,2\sqrt{P}$ $d = 1,2\sqrt{P} \text{ 80 MHz thru 800MHz}$ $d = 2,3\sqrt{P} \text{ 800 MHz thru 2,5MHz}$ Where P is the nominal maximum power of output of the transmitter in watts (W), as per the manufacturer of the transmitter, and d is the recommended separation distance in meters (m). It is advisable that the field intensity from the RF, transmitter as determined by means of electric inspection on-site, ^a is less than the level of compliance in each frequency range ^b. There may be interference near the equipment marked with the following symbol: 
NOTE 1 At 80MHz and 800MHz, the highest frequency range applies.			
NOTE 2 These directives may not be applicable in every situation. The electromagnetic transmission is affected by the absorption and reflection of structures, objects and people.			
^a The field intensities set by the fixed transmitters, such as radio base stations, telephones (mobile phone, wireless) land mobile radio, amateur radio, AM and FM radio transmissions and TV transmissions can not be predicted with accuracy. Due to the RF fixed transmitters is recommended to install an electromagnetic inspection at the local in order to evaluate the electromagnetic environment. If at the place where the equipment is be using the field intensity level exceeds the conformity level for the RF above, is recommended to observe if the operations are normal. Whether abnormal operations are observed, additional procedures shall be necessary such as reorientation or replace the equipment. ^b Whether above the frequency range of 150kHz to 80 MHz is recommended a field intensity below than 3 V/m.			

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TECHNICAL SPECIFICATIONS

Eletromagnetic emissions		
The Celeritá Electrical Micromotor is made to be used in the electromagnetic environments specified below. The client or the user of the Celeritá Electrical Micromotor must be sure that it is used in such environment.		
Emission test	Compliance	Eletromagnetic environment - Guide
RF emissions ABNT NBR IEC CISPR 11	Group 1	This equipment uses RF energy only for internal functions. However, its emissions are too low and it's unlikely to cause any interference in the equipments next to it.
RF emissions ABNT NBR IEC CISPR 11	Class B	This equipment is proper to be used in all establishments; including domestic settings and those directly connect to a public low voltage distribution which feeds domestic buildings.
Emissions of harmonics IEC 61000-3-2	Class A	
Fluctuation of Voltage / Emissions of flicker IEC 61000-3-3	As per	

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TECHNICAL SPECIFICATIONS

The materials used to produce the equipment are Biocompatible


The equipment pieces are not subject to degradation, corrosion, fatigue or material aging, during its shelf life.







List of pieces and circuit scheme

Gnatus Company declares that the supply of the circuit scheme, list of pieces or any other information that propitiates technical attendance for the user, can be requested if there is an agreement between the user and Gnatus Company.



Dimensions (mm)

Packing symbols

Warning - see the manual

Content of accessible and non-accessible demarcations

OPERATION OF EQUIPMENT

General use

A) To fix the drill, ensure that the pliers is open. Rotate the regulator ring clockwise (position R) until you hear a click.

B) Insert the drill completely; after touching the bottom of the tweezer, return the drill for about 2 mm. Rotate the regulator ring in the contrary direction (position S) blocking the drill.

! For your safety, Micromotor has factory default with the voltage selector switch set to 220V, so make sure that your network is compatible before turning it on; if it is on 127V, reverse the position.

ATTENTION:

Do not reverse the voltage with the device turned on.

- Place the key (C) according to the Local voltage (127 or 220V).

Afterwards turn the equipment on and proceed according to the sequence of operations as follows:

- Connect the hand piece cable to "MOTOR" connection (E) of the control unit;
- Turn the switch on (F);
- Select in the switch (G) the option hand (H);
- Select in the switch (D) the position "FWD"(forward) to turn the engine in only one direction, or "REV" to the reverse direction;

The engine speed can vary between 0 and 20.000 rpm through the controller button of the potency (I).

OPERATION OF EQUIPMENT

IMPORTANT:

- When the handpiece is in operation (turning), do not try to adjust the pliers, in order to prevent damage to the control unit;
- Before operating the FWD/REV button, be sure that the engine is turned off;
- Do not leave the equipment turn on with the potentiometer in low potency (Low), in order to prevent engine overheating.

PRECAUTIONS, RESTRICTIONS AND WARNINGS

Transportation, storage and operation

This equipment must be transported and stored observing the following directions:

- Avoid falls and impacts;
- Keep it dry, do not expose it to rain, water drops or wet floor;
- Keep it away from water and direct sunlight, and in its original wrapping;
- Don't move it over irregular surfaces, protect it from rain and observe the maximum stack quantity specified in the packaging;
- Transportation and storage temperature range: -12°C to 50°C.

The Equipment maintains its condition of safety and efficacy, provided that it is maintained (stored) as mentioned in this instruction of use. Thus, the equipment will not lose or alter its physical and dimensional features.

Sensitivity to environmental conditions in normal situations of use

- The equipment has been planned not to be sensitive to interference such as magnetic fields, external electrical factors, electrostatic discharge, pressure or variance of pressure, provided that the equipment is installed, maintained, clean, preserved, transported and operated as per this instruction for use.

Precautions and warnings “during the installation” of equipment

- Check the voltage of the equipment at the moment of executing the electrical installation.
- Position the equipment in a place where it will not get wet.
- Install the equipment in a place where it will not be damaged by pressure, temperature, humidity, direct sunlight, dust, salts, or sulfur compounds.
- The unit should not be submitted to inclination, excessive vibrations, or blows (including during transportation and handling).
- This equipment was not planned for use in an environment where vapors, anesthetic mixtures inflammable with air, or oxygen and nitrous oxide can be detected.
- Before the first use and/or after long interruptions from work such as vacations, clean and disinfect the equipment..

PRECAUTIONS, RESTRICTIONS AND WARNINGS

Recommendations for the dental equipment maintenance

Your Gnatus equipment has been designed and developed according to the standards of modern technology. Similarly to other kinds of equipment, it requires special care, which is many times neglected due to several reasons and circumstances.

Therefore, here are some important reminders for your daily routine. Try to follow these simple rules, which will save you a lot of time and will avoid unnecessary expenses once they start making part of your working procedure.

Precautions and warnings “during the use” of equipment

- The device shall be operated only by qualified personnel properly authorized and trained (Qualified Professionals)

- In the need of maintenance, use only the services of Gnatus Authorized Technical Assistance.

- Although this device has been designed according to electromagnetic compatibility standards it may, under extreme conditions, cause interference with other devices. Do not use this device with other devices that are very sensitive to interference or with devices that create high electromagnetic disturbance.

- Never block the openings of ventilation of the device. Keep the openings free from Lints, hairs and other materials.

- Never turn the equipment on without the drill. It may cause pliers blocking.

- Be always alert to the quality and wearing of drill. Drills in good conditions reduce the work time resulting in less wear of bearings.

- After replacing the drill, it is necessary ensure that it is fixed to the pliers.

- If the drill is used while is not fully inserted, it is dangerous and may be loose due to its vibration.

- Do not turn the ring that releases the pliers of Micromotor while it is in operation.

- Never use twisted, asymmetric or damaged drills.

Gnatus shall not be responsible for:

- Use of the equipment differing from that for which it is intended.

- Damages caused to the equipment, the professional and/or the patient by the incorrect installation and erroneous procedures of maintenance, differing from those described in these Instructions for use which come with the equipment or by the incorrect operation of it.

Precautions and warnings “after” the use of equipment

- Turn off the main switch of the dental set when it is not in use for an extended period of time.

- Always maintain the equipment clean for the next operation.

- Do not modify any part of the equipment. Do not disconnect the cable or other connections without need.

- After using the equipment, clean and disinfect all the parts which may be in contact with the patient.

- Always check the unit in relation to its looseness, vibration, sound or temperature in the operation. If some abnormal operation is detected, automatically stop the operation and contact Gnatus Technical service.

PRECAUTIONS, RESTRICTIONS AND WARNINGS

Precautions and warnings during the “cleaning and disinfection” of equipment

- Before cleaning the equipment, turn off the main switch.
- The hand piece cannot be washed, sterilized, autoclaved or lubricated because Shielded bearings have special grease and may be damaged.
- Avoid spilling water, even accidentally, or other liquids inside the equipment, which could cause short circuits.
- Do not use microabrasive material or steel wool when cleaning, or employ organic solvents or detergents which contain solvents such as ether, stain remover, gasoline etc.

Precautions in case of alteration in the functioning of equipment

- If the equipment has any abnormality, check if the problem is related to any item listed in the topic of unforeseen events (failures, causes and solutions). If it is not possible to resolve the problem, turn off the equipment, remove the power supply cable from the socket and contact your representative (Gnatus).

Precautions to be adopted against foreseeable or uncommon risks, related to the deactivation and abandoning of equipment

In order to avoid environmental contamination or undue use of the Equipment after it has become useless, it should be discarded in the suitable place (as per the local legislation of the country).

- Pay attention to the local legislation of the country for the conditions of installation and disposal of residue.

CORRECTIVE AND PREVENTIVE MAINTENANCE AND PRESERVATION

Preventive Maintenance

In order to protect your equipment, seek Gnatus technical assistance for periodic revisions of preventive maintenance.

Corrective Maintenance

If the equipment has any abnormality, check if the problem is related to any of the items listed in the item Unforeseen Events (situation, cause and solution).

If it is not possible to solve the problem, turn off the equipment, and request Gnatus technical assistance.

Additional procedures for reuse

The equipment can be reused in undetermined, i.e. unlimited, quantities, only needing to be cleaned and disinfected.

CORRECTIVE AND PREVENTIVE MAINTENANCE AND PRESERVATION

Cleaning

Important: In order to execute cleaning or any type of maintenance, ensure that the equipment is disconnected from the electrical network.

The cleaning procedure should be executed at the start of the working day and after each patient. Always turn the main switch off before executing the procedures of daily maintenance.

To clean the equipment, we recommend the use of "BactSpray (Reg n° MS: 3.2079.0041.001-5) or any other similar product:

Active component: Benzalkonium chloride (tri-quaternary ammonium)


Solution 50%..... 0.329%

Chemical composition: Butyl Glycol, Decyl polyglucose, Sodium Benzoate, Sodium Nitrate, Essence, Deodorized Propane / Butane, demineralized Water.

For more information concerning cleaning procedures, see manufacturer's instructions.

WARNING:

- In order to prevent risks and damages to equipment, make sure that the liquid does not enter into the unit.
- The application of other solvent-based cleaning products or sodium hypochloride isn't recommended, because they may damage the equipment.

NOTE: The registration at the Ministry of Health of the "BactSpray" is executed separately from the product described in this manual, as the "BactSpray" is not manufactured by Gnatus.

NOTE: Use gloves and other systems of protection, during the disinfection.

Disinfection

- Use clean and soft cloth dampened in alcohol 70% to disinfection of the equipment.
- Never use corrosive disinfectants or solvents.

UNFORESEEN EVENTS – SOLUTION OF PROBLEMS

 Upon coming across any problem in operation, follow the instructions below to check and repair the problem, and/or get in touch with your representative.

Problem	Probable cause	Solution
-Inoperable micromotor.	-Main Switch in turn off position. -Disconnected Plug from network. -No power. -Failure in the connection between micromotor cable and control unit. - Damaged Items: -Micromotor cable; -Fuse; -Micromotor.	-Turn the main switch on. -Connect the plug. -Wait for network normalization. -Connect correctly the cable of the Micromotor. -Contact Gnatus technical service.
-Micromotor suddenly stops during use. "Protection system operated."	-Overcharge. -Open pliers during operation.	-See item: "safety and protection sy - pag.23". -See the pliers and close it.
-Micromotor is not operating when the pedal is activate.	-Selector hand/foot is in "H" (hand). -Foot control operating failure. -Foot control cable damaged. -Failure connection between the pedal cable and control unit.	-Place the selector in the "F" position (foot). -Contact Gnatus technical service. -Contact Gnatus technical service. -Check connection.

WARRANTY OF EQUIPMENT

This equipment is covered by the warranty terms counting from the date of installation, as specified below; provided that the defect has occurred in normal conditions of use and that the equipment has not remained stored for more than 06 months counting from the issue date of the sales document until the date of the actual installation.

- WARRANTY TERMS: TWELVE MONTHS;

- LOSS OF THE WARRANTY:

- A) Attempt to repair using an inadequate tool or by unauthorized technicians;

- B) Installation of the equipment by an unauthorized technician;

- C) Damage arising from inappropriate storage or signs of infringement;

- D) Incorrect use of the equipment;

- E) Use of a cleaning product not indicated by the factory;

- F) Falls or blows which the equipment may undergo or lack of observation of an compliance

with the guidelines of the Owner's Manual, which was delivered with the present document, together with the equipment. Repair or replacement of parts during the warranty period shall not extend the validity term of their warranty.

- This warranty does not exempt the customer from paying the service charge for the visit and the travel expenses of the technician, except when the customer sends the equipment to execute the maintenance inside the establishment of the technical assistance.

- "Consumer Defense Code - art. 50, unique paragraph".

- The Warranty Certificate comes with the product and must be filled in upon the date

- of installation by the Gnatus Authorized Technician.

- Queries and information: GNATUS Help Desk (+55) 16 2102-5000.

FINAL CONSIDERATIONS

The most important aspect related to equipment care is that concerning spare parts.

To guarantee the life span of your equipment, use only original Gnatus spare parts. They are sure to follow the technical specifications and standards required by Gnatus.

We must also point out to you our chain of authorized dealers. Only dealers that make part of this chain will be able to keep your equipment constantly new for they count on technical assistants who have been trained and on specific tools for the correct maintenance of your equipment.

Doubts and information: GNATUS Call center (55-16) 2102-5000.

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